

B'SYS GmbH

Kainate-Receptor Assay

Overview Sheet

© B'SYS GmbH

TABLE OF CONTENTS

3.	.1. Contact Address for Technical Support & Ordering Information	5
3.	CONTACT INFORMATION	5
2.	ZK200775, A COMPTITIVE ANTAGONIST	4
1.	GLUTAMATE DEPENDENCE	3



1. GLUTAMATE DEPENDENCE

Briefly, CHO cells were transfected with GluK2 or GluK2/5 and split for electrophysiological experiments after 24-48 h. During the electrophysiological experiment they were continuously perfused with buffer solution: 137 mM NaCl, 4 mM KCl, 1.8 mM CaCl $_2$, 1 mM MgCl $_2$, 10 mM Glucose, 10 mM HEPES. pH 7.4 (NaOH). Pipette was filled with a solution containing: 130 mM KCl, 1 mM MgCl $_2$, 5 mM EGTA, 10 mM HEPES, pH 7.2 (KOH). After the seal was obtained cells were clamped to -80 mV. GluK2 channels were activated with increasing concentrations of Glutamate 10 μ M - 3000 μ M (application interval: 30 s).




Figure 1: Effect of Glutamate on (A) GluK2 homomeric and (B) heteromeric GluK2/5 channels.

Table 1: Glutamate Dependence.

Receptor composition	EC ₅₀	Hill coefficient
GluK2 homomer	159.32 ± 7.32 μM	3.58
GluK2 / GluK5 heteromer	3.47 ± 0.63 μM	1.37



2. ZK200775, A COMPTITIVE ANTAGONIST

ZK200775 is a competitive antagonist for that displays anxiolytic, anticonvulsant and muscle relaxant activity *in vivo*. Concentrations between 1.0 and 400 μ M were tested in the presence of 80 μ M Glutamate for GluK2 or 5.0 μ M for GluK2/5.




Figure 2: Effect of ZK200775 on (A) GluK2 homomeric and (B) heteromeric GluK2/5 channels. Application scheme: Glutamate, Glutamate + 1.0 μ M ZK200775, Glutamate + 10 μ M ZK200775, Glutamate + 400 μ M ZK200775

Table 2:

Receptor composition	IC ₅₀	Hill coefficient
GluK2 homomer	38.24 μM	0.66
GluK2 / GluK5 heteromer	21.12 μM	0.88



3. CONTACT INFORMATION

3.1. Contact Address for Technical Support & Ordering Information

B'SYS GmbH
 Technology Center Witterswil
 Benkenstrasse 254
 4108 Witterswil
 Switzerland

Tel: +41 61 721 77 44
Fax: +41 61 721 77 41
Email: assay@bsys.ch
Web: www.bsys.ch

